

Extra Heavy Duty Body Bag with Curved Zipper

Extra Heavy Duty Body Bag with Curved Zipper

Designed by Basic Medical Supply. Body bag has 6 strong handles, 3 each side. This is commonly referred to as a Human Remains Pouch, Human Remains Bag, Human Remains Pouch or a Body Bag. Material: 12 mil PVC. Static Lift Tested to 300 LB Capacity. Each Bag contains 3 Toe Tags. All Packed 10/case.

ITEM #	DESCRIPTION	QTY	PRICE
17014	36" x 96" Color: Black Material: 20Mil PVC	Each	

