

**Basic Medical Supply ( CHINA)** 11-1703, # 225 Jinxiang Road Pudong, Shanghai, China Tel: +86 (21) 3175 0988 Fax: +86 (21) 3175 0989

**Basic Medical Supply ( USA )** 19902 Flax Flower Dr. Richmond, TX 77407 Tel:+1-832-535-3107 Fax:++1-832-535-3107

www.basic-medicalsupply.com

## **Clear Tape**

## **Clear Tape**

Hypoallergenic tape is transparent to allow for visualization of I.V. tubing, catheters, and other dressings that require monitoring and features perforations to aid in air circulation. This tape tears easily both across and lengthwise.

ITEM #	DESCRIPTION	QTY	PRICE
MS-TT1/2	1/2" x 10yds	Each	
MS-TT1	1" x 10yds	Each	
MS-TT2	2" x 10yds	Each	
MS-TT3	3" x 10yds	Each	
MS-TT1225	1/2" x 2.5yds	Each	
MS-TT15	1" x 5yds	Each	



info@basic-medicalsupply.com

www.basic-medicalsupply.com

