

## Deluxe Evacuation Stair Chair(401)

### Deluxe Evacuation Stair Chair(401)

The Deluxe Evacuation Stair Chair is tough and lightweight and folds to a depth of only 8" for easy vehicle transport. The chair features extended lift handles on the foot end for safer lifting, better line of sight and increased distance between patient and operator. The chair is powder-coated and features front caster wheels and the stair-tread system to control movement in stairways without lifting. The chair is an essential tool for evacuating disabled or injured patients from multi-level buildings. The chair measures 37.5" x 20.5" x 28" and weighs 31.5 lbs.

The folded depth is 8" and the maximum weight capacity is 500 lbs.

ITEM #	DESCRIPTION	QTY	PRICE
MS-STA411	Yellow, Aluminum	Each	

