

Folding Stretcher, Aluminum

Folding Stretcher, Aluminum

Folding stretcher strengthened by two steel bars that allow its use even for heavy transports. It is foldaway for an easy location. Orange, high quality PVC sheet is characterized for its fire proof, and water proof, anti-tear and easy to clean. Product

Size (L*W*H):2060*520*140mm

Folded Size (L*W*H):1030*180*100mm

Carton Size (2pcs): 1130*210*210mm

Load Limit: 159kg

N.W.:5kgG.W.: 11kg

ITEM #	DESCRIPTION	QTY	PRICE
MS-FL002	Orange, Size: 2080*545*130mm, Material: Aluminum	Each	

