

Basic Medical Supply ( CHINA) 11-1703, # 225 Jinxiang Road Pudong, Shanghai, China Tel: +86 (21) 3175 0988 Fax: +86 (21) 3175 0989 **Basic Medical Supply ( USA )** 19902 Flax Flower Dr. Richmond, TX 77407 Tel:+1-832-535-3107 Fax:++1-832-535-3107

www.basic-medicalsupply.com

## **Folding Stretcher, Aluminum**

## Folding Stretcher, Aluminum

Folding stretcher strengthened by two steel bars that allow its use even for heavy transports. It is foldaway for an easy location. This stretcher is intensified stretcher, with flexible handle and outrigger, can be folded by lengthwise. It characterized by its being easy-carried, use-safely and easy for sterilization and clean. Product Size (L\*W\*H):2100\*550\*120mm Folded Size (L\*W\*H):1020\*550\*160mm Carton Size (1pc): 103\*550\*180mm Load Limit: 159kg N.W.:10kg G.W.: 11kg

ITEM #	DESCRIPTION	QTY	PRICE
MS-FL009	Orange, Size:2100*550*120mm, Material: Aliuminum	Each	



info@basic-medicalsupply.com

**basic** MEDICAL SUPPLY

Basic Medical Supply (CHINA) 11-1703, # 225 Jinxiang Road Pudong, Shanghai, China Tel: +86 (21) 3175 0988 Fax: +86 (21) 3175 0989 Basic Medical Supply ( USA ) 19902 Flax Flower Dr. Richmond, TX 77407 Tel:+1-832-535-3107 Fax:++1-832-535-3107

www.basic-medicalsupply.com