

Basket Stretcher(2000)

Basket Stretcher(2000)

The stretcher was made with an eye on a very wide range of emergency situations and ,as it is sturdy and flexible to use,it is a must for all types of aid in special situations,it's quick reliable fittings enable first-aid personnel to operate rapidly and safely. With its special sling equipment,the stretcher is ideal for lifting and transport by helicopter,it has an adjustable feet-securing mechanism,safety belts and a mattress.It is also designed for to be secured inside the helicopter,All the materials used have dimensions suitable to work and they are fire-proof,they do not release any toxic

or polluting substances and they are protected to obtain high resistance to wear and corrosion. Product Size (L*W*H): 2190*640*180mm

Carton Size: 1300*680*310mm

Load Limit: 272kg N.W.:11kg G.W.:20kg

Item # MS-BS2000

ITEM #	DESCRIPTION	QTY	PRICE
MS-BS2000	Orange, Size: 2190*640*180mm	Each	

