

6 Handles Sit Stretcher

6 Handles Sit Stretcher

The Sit Stretcher is designed to transport patients from conined areas. Patient is carried in an upright sitting position. Stretcher is compact measuring 20" x 34" and is constructed of easy to clean

orange vinyl and 2" continuous yellow nylon straps with 6 handholds. Easy to clean orange vinyl carrying case is included.

ITEM #	DESCRIPTION	QTY	PRICE
22017BL	Orange; 20"x34"; Impervious Fabric; 20pcs/Case	Each	
22017OR	Orange; 20"x34"; Impervious Fabric; 20pcs/Case	Each	

